

A decontamination solution created specifically for Getz

Supporting Getz with their healthcare systems.



Providing valuable healthcare solutions to individuals in the Asia Pacific region.

A leading distributor of medical devices, equipment, and consumables in Hong Kong since 1978.

Getz Healthcare (Hong Kong) Ltd operates in various healthcare fields, such as Interventional Cardiology, Radiology, Infection Control, Pulmonology, Minimally Invasive Surgery, and Vascular Surgery.

Getz has been providing professional services and experienced personnel across Hong Kong SAR, Macau, and other Asian countries since 1978, representing 24 brands.

Covering many areas of medical practice with quality products from world-class manufacturers, serving their customers with first-class and value-added services within:

- Bioanalytics
- Dental
- Diagnostics
- Medical

PBSC provided high-quality decontamination solutions using chlorine dioxide (ClO₂) and updated software compatible with ClO₂ generators, supporting Getz and their healthcare systems.

Products


- PBSC - 2 Valve Air Handling System design
- H14 HEPA Filtration for Chamber inlet air
- Pressure relief valve
- Chlorine Dioxide Generator
- Gas Tight Damper Valves



Solutions/Results

PBSC provided a decontamination solution that can accommodate Getz's preferred method of using Chlorine Dioxide, with high quality and flexibility.

- **PBSC's 2 Valve Air Handling System**
 - Allowing to duct their HVAC system straight from/to PBSC's inlet and extract valves.
- **Gas Tight Damper Valves**
 - The inlet and outlet valves were ordered with the gas tight option, which offers a higher level of containment.
- **Remote control panel and pneumatics panel**
 - Getz is able to install these products into technical areas more easily due to their remote capabilities, away from production and higher classification areas.



- **HVAC and BMS I/Os**
 - PBSC offered volt-free signals, available if Getz requires the notification of alarms via their BMS system to indicate status.
- **Pressure Relief Valve**
 - This valve, in conjunction with the pressure relief HEPA/CAT filter pack, enables pressure testing and bleeding of pressure during a cycle.
- **Chlorine Dioxide Generator**
 - The PBSC Chamber allows Getz to operate their Chlorine Dioxide Generator from the Chamber HMI screens via volt free contacts, making it highly versatile.

PBSC provided essential decontamination products to aid Getz in their comprehensive service and solutions.



“ Emma Gibbons
- Senior Project
Manager for PBSC
- We are pleased
to have provided
Getz with our
decontamination
solutions and wish
them every success
in their future
endeavors ”

www.pbsc.co.uk
www.pbsc-inc.com